

Noise Equivalent Power (NEP)

Electronic circuits output noise as well as the desired signal; an example would be the hiss heard in a stereo when no music is playing. The noise spectrum has a relatively flat response, and the noise level changes with the square root of the frequency range; for example, if the frequency range doubles, the noise component increases by $\sqrt{2}$ (1.414). The optical input has to be large enough to overcome the noise component of the detector in order to have a measurable signal output. Noise equivalent power (NEP) is a specification that allows a customer to determine the noise component of a detector for their particular wavelength and measurement bandwidth. NEP is the minimum optical power required for an output signal-to-noise ratio (SNR) of 1. This means the signal level and the noise level are the same; an SNR of 10 or larger is recommended. NEP has units of $W/\sqrt{\text{Hz}}$ or W , depending on whether or not the measurement bandwidth is included.

NEP is calculated from measurement or photodiode parameters such as dark current (the current that flows through a reverse biased photodiode when light is not present). Typical values and noise sources are:

Non-amplified standard detector: $0.1\text{pW}/\sqrt{\text{Hz}}$ (dark current)

High speed detector: $30\text{pW}/\sqrt{\text{Hz}}$ (noise from internal 50Ω impedance matching resistor)

Amplified detector: $25\text{pW}/\sqrt{\text{Hz}}$ (amplifier noise)

EOT's data sheets provide NEP at a specific wavelength – this value is needed to determine NEP for other wavelengths. An example is shown below for the ET-2030:

Measurement Wavelength: 532nm

Measurement Bandwidth: 800MHz

Responsivity (R_λ) of ET-2030 at 532nm (from graph in data sheet): 0.27A/W

Responsivity (R_λ) of ET-2030 at 830nm (from table in data sheet): 0.47A/W

Datasheet NEP: $0.01\text{pW}/\sqrt{\text{Hz}}$

$$\text{NEP}_{532\text{nm}} = \text{NEP}_{\text{datasheet}} \times \frac{R_{\text{datasheet}}}{R_{532\text{nm}}} = \frac{0.01\text{pW}}{\sqrt{\text{Hz}}} \times \frac{0.47\text{A/W}}{0.27\text{A/W}} = 0.02\text{pW}/\sqrt{\text{Hz}}$$

Including the bandwidth:

$$\text{NEP}_{532\text{nm}} = \frac{0.02\text{pW}}{\sqrt{\text{Hz}}} \times \sqrt{(800 \times 10^6\text{Hz})} = 566\text{pW}$$

For an SNR of 10 the optical power would need to be 5.66nW (-52.5dBm). This value would be fine for a spectrum analyzer, but oscilloscopes typically require an output of 1mV or higher for a minimum value. In practice, this is too low of a value: 10mV to 20mV would be an acceptable



minimum reading. Using the ET-2030 as an example, and working backward from the oscilloscope output of 1mV:

$$\text{Optical Input} = \frac{\text{Output Voltage}}{\text{Scope Termination}} \times \frac{1}{R_{532\text{nm}}} = \frac{0.001\text{V}}{50\Omega} \times \frac{\text{W}}{0.27\text{A}} = 74\mu\text{W}$$